

Product Name Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 5 (ST6GALNAC5)

Catalog Number #0026

Alternate Names GD1 alpha synthase, GalNAc alpha-2,6-sialyltransferase V, ST6GalNAc V, (ST6GalNAcV),

Sialyltransferase 7E (SIAT7-E)

Substrate Specificity redominantly catalyzes the biosynthesis of ganglioside GD1alpha from GM1b in the brain, by

transferring the sialyl group (N-acetyl-alpha-neuraminyl or NeuAc) from CMP-NeuAc to the

GalNAc residue on the NeuAc-alpha-2,3-Gal-beta-1,3-GalNAc sequence of GM1b Furukawa et al., 'ST6 N-acetylgalactosaminide Alpha-2,6- sialyltransferase 5,6

(GST6GALNAC5,6)' in Handbook of Glycosyltransferases and Related Genes, 2nd Edition

Expression Host HEK293
Species of expressed protein Human
Gene ID 81849
Protein RefSeq NP 112227

 Uniprot
 Q9BVH7

 Region Expressed
 AA 50 - 336

Expressed Protein Sequence

References

ASATGSSQPAAESSTQQRPGVPAGPRPLDGYLGVADHKPLKMHCRDCALVTSSGHLLHSRQ GSQIDQTECVIRMNDAPTRGYGRDVGNRTSLRVIAHSSIQRILRNRHDLLNVSQGTVFIFWGPS SYMRRDGKGQVYNNLHLLSQVLPRLKAFMITRHKMLQFDELFKQETGKDRKISNTWLSTGWFT MTIALELCDRINVYGMVPPDFCRDPNHPSVPYHYYEPFGPDECTMYLSHERGRKGSHHRFITE

KRVFKNWARTFNIHFFQPDWKPESLAINHPENKPVF

Tag(s) N-terminal 6xHis, GFP


Purity (%) >95%, by SDS-PAGE as visualized by Coomassie Blue Staining

Formulation Supplied as a 0.2µm filtered solution in 20mM HEPES pH 7.0, and 100mM NaCl buffer, with

10% Glycerol and 0.05 % NaN₃ as preservative.

Concentration 1 mg/ml SDS-Page Size ~70 kDa

SDS-PAGE image



Activity Measured by the ability to transfer the sugar from CMP-Neu5Ac and generate CMP

Assay Buffer 100mM MES, pH 7.0

Donor Substrate CMP-Neu5Ac (200 mM, Nacalai Tesque Inc.)

Acceptor Substate L

Detection Kit Assay Steps LNT

CMP-Glo™ Glycosyltransferase Assay (Promega)

- 1) Prepare 10ml reaction mixture containing 100mM MES (pH 7.0), CMP-Neu5Ac (200mM) as donor and LNT (1 mM) as acceptor and purified GFP-ST6GALNAC5 in a microfuge tube.
- 2) Incubate at 37C° for 30 min.
- 3) Put the sample on ice immediately and then transfer 5 µLof reaction mixture into 384-well assay plates and add equal volume of CMP Detection Reagent (5µL)

4) Incubate for 60 min at room temperature and read the plate using a GloMax Multi Detection System plate reader (Promega)

Std Curve Follow protocol for "Generating a Standard Curve for CMP" in the CMP-Glo™

Glycosyltransferase Assay Technical Manual (Promega)

Specific Actifity calc Specific Activity (pmol/min/ug)= [CMP released*(nmol) x (1000 pmol/nmol)] / [Incubation time

(min) x amount of enzyme (ug)], Specific Activity was calculated using the standard curve

plotted in GraphPad Prism 6 (GraphPad Software)

Shipping conditions This product is shipped as 0.2µm filtered product on dry ice. Upon receipt, store it immediately

at the temperature recommended below.

Stability & Storage cond 6 months if stored at -80C. Avoid repeated freeze thaws.